

Product Data Sheet

Indoxyl-Glucoside

Product No. GB-CHEM-00625
Volume. 100MG

Quality Profile	Specification
Appearance (Color)	White to Off-White
Appearance (Form)	Powder
Solubility (Color)	Colorless to Yellow
Solubility (Turbidity) 50 mg/mL in H ₂ O	Clear
Water By Karl Fischer	≤ 10
Specific Rotation (anhydrous) At 546 nm, C = 1 in H ₂ O at 20 °C	-80 - -73 °
Purity (TLC)	≥ 97 %

Intended Use: For research use only.

This document has been created digitally and is valid without a signature.