

Product Data Sheet

Kinetin

Product No. GB-CHEM-00626
Volume. 100G

Quality Profile	Specification
Appearance (Color)	White to Light Beige
Appearance (Form)	Powder to crystals
Purity (HPLC area %)	≥ 99.0 %
Melting point	269 - 271 °C
Residue on ignition	≤ 0.05 %
Carbon content	55.81 % (theory)
Hydrogen content	4.22 % (theory)
Nitrogen content	32.54 % (theory)
Infrared spectrum	Conforms to Structure

Intended Use: For research use only.

This document has been created digitally and is valid without a signature.