



Product Data Sheet

Lysostaphin

Product No. GB-CHEM-00629
Volume. 1MG

Quality Profile	Specification
Appearance (Color)	Lyophilized Powder
Appearance (Form)	White to Light Yellow
Solubility (Solvent)	Water
Solubility (Conc)	2.00 - 2.10 mg/mL
Solubility (Turbidity)	Clear to Slightly Hazy
Solubility (Color)	Colorless to Light Yellow
Protein Content (Method)	Biuret
Prot. Content (% Protein)	50.00 - 70.00 %
Enzym. Assay (Activity U/mgS)	≥ 0.00
Enzym. Assay (Spec. Activity)	≥ 500.00

Intended Use: For research use only.

This document has been created digitally and is valid without a signature.