



Product Data Sheet

Malachite green oxalate

Product No. GB-CHEM-00631
Volume. 50G

Quality Profile	Specification
Appearance (Color)	Faint Green to Very Dark Green
Appearance (Form)	Powder or Solid or Crystals or Chunk(s) Crystalline Chunk(s) or Chunk(s)
Solubility (Color)	Blue to Very Dark Blue
Solubility (Turbidity) 1 mg/ml, H ₂ O	Solution
Infrared Spectrum	Conforms to Structure
Dye Content	≥ 90 %
Biological Stain Commission	Pass

Intended Use: For research use only.

This document has been created digitally and is valid without a signature.