



## Product Data Sheet

### Manganese sulfate, monohydrate

**Product No.** GB-CHEM-00632  
**Volume.** 500G

Quality Profile	Specification
Appearance (Color)	Faint Red
Appearance (Form)	Powder
Solubility (Color)	Almost Colorless to Faint Red
Solubility (Turbidity) 10 g/130 mL, H <sub>2</sub> O	Clear
Identity - Manganese (Mn)	Pass
Titration with EDTA	≥ 99 %

**Intended Use:** For research use only.

This document has been created digitally and is valid without a signature.