



Product Data Sheet

Ninhydrin

Product No. GB-CHEM-00635
Volume. 100G

Quality Profile	Specification
Appearance (Color)	White to Yellow and Faint Beige to Light Beige and Faint Brown to Light Brown
Appearance (Form)	Powder
Suitability Clear dark yellow to green Solution at 400 mg plus 4 ml of Ethanol	Pass
Loss on Drying	≤ 13.2 %
Carbon Content	59.5 - 61.9 %

Intended Use: For research use only.

This document has been created digitally and is valid without a signature.