



Product Data Sheet

Nystatin

Product No. GB-CHEM-00636
Volume. 50G

Quality Profile	Specification
Appearance (Color)	Yellow
Appearance (Form)	Powder
Suitability	Corresponds
Water (by LOD)	≤ 6.0 %
pH 3% solution	6.0 - 8.0
Potency (On Dry Basis) USP unit/mg	≥ 4400
Infrared spectrum	Conforms to Structure
Cell Culture Test	Pass
Plant Cell Culture Test	Pass

Intended Use: For research use only.

This document has been created digitally and is valid without a signature.