



Product Data Sheet

Rubidium chloride

Product No. GB-CHEM-00641
Volume. 25G

Quality Profile	Specification
Appearance (Color)	White to Off White
Appearance (Form)	Powder or Crystals
Solubility (Colour)	Colorless to Very Faint Yellow
Solubility (Turbidity) 100 mg/mL in water	Clear to Slightly Hazy
UV Absorbance 280nm (50% W/V IN WATER)	≤ 0.4
Metal Trace Analysis (ICP)	≤ 10000 ppm
Purity	≥ 99.0 %

Intended Use: For research use only.

This document has been created digitally and is valid without a signature.