



Product Data Sheet

Salicylic acid

Product No. GB-CHEM-00642
Volume. 250G

Quality Profile	Specification
Appearance (Color)	White
Appearance (Form)	Powder or Crystals
Melting Point	158.0 - 161.0 °C
Infrared Spectrum	158.0 - 161.0 °C
Purity (HPLC)	≥ 99.0 %
Residue on ignition (Ash)	≤ 0.01 %
Chloride Content	≤ 0.001 %
Substances Darkened (by H ₂ SO ₄)	Pass
Iron (Fe)	≤ 2 ppm
Heavy Metals (as Lead)	≤ 5 ppm
Sulfate (SO ₄)	≤ 0.003 %

Intended Use: For research use only.

This document has been created digitally and is valid without a signature.