



Product Data Sheet

ACES

Product No. GB-CHEM-00646
Volume. 25G

Quality Profile	Specification
Appearance (Color)	White
Appearance (Form)	Crystalline Powder
Solubility (Color)	Colorless
Solubility (Turbidity) 5 g plus 95 ml H ₂ O	Clear
Infrared Spectrum	Conforms to Structure
Loss on Drying	≤ 1.0 %
A ₂₉₀ UV absorbance (c = 5% w/w)	≤ 0.02
Titration with NaOH (dry basis)	≥ 99.0 %

Intended Use: For research use only.

This document has been created digitally and is valid without a signature.