



Product Data Sheet

Acetic Acid (Glacial)

Product No. GB-CHEM-00648
Volume. 1L

Quality Profile	Specification
Assay (acidimetric)	100.0%
Color	≤ 10 Hazen
Titrateable base	≤ 0.0004 meq/g
Acetic anhydride	≤ 100 ppm
Chloride (Cl)	≤ 1 ppm
Heavy metals (as Pb)	≤ 0.5 ppm
Sulfate (SO ₄)	≤ 1 ppm
Fe (Iron)	≤ 0.2 ppm
Evaporation residue	≤ 10 ppm

Intended Use: For research use only.

This document has been created digitally and is valid without a signature.