

Product Data Sheet

Acetone

Product No. GB-CHEM-00649
Volume. 1L

Quality Profile	Specification
Appearance (Color)	Colorless
Appearance (Form)	Liquid
Infrared Spectrum	Conforms to Structure
Corrected for Water Content Purity by GC	≥ 99.50 %
Methanol Content	≤ 0.05 %
GC Identified Peak Isopropyl Alcohol	≤ 0.05 %
Water (by Karl Fischer)	≤ 0.5 %
Residue on Evaporation	≤ 0.001 %
Color Test	≤ 10 APHA
Solubility	Pass
Aldehyde (as Formaldehyde)	≤ 0.002%
Substances reducing MnO ₄	Pass
Titrateable Acid (meq/g)	≤ 0.0003
Titrateable Base (meq/g)	≤ 0.0006

Intended Use: For research use only.

This document has been created digitally and is valid without a signature.