



Product Data Sheet

3',5'-Dimethoxy-4'-hydroxyacetophenone

Product No. GB-CHEM-00651
Volume. 1G

Quality Profile	Specification
Appearance (Color)	White to Yellow and Faint Beige to Beige Off White and Faint Brown to Brown
Appearance (Form)	Powder or Crystal or Chunk(s)
Titration with TBAOH	96.5 - 103.5 %
Purity (HPLC)	≥ 96.5 %
Infrared Spectrum	Conforms to Structure

Intended Use: For research use only.

This document has been created digitally and is valid without a signature.