



## Product Data Sheet

### **Acriflavine, neutral**

**Product No.** GB-CHEM-00656  
**Volume.** 25G

<b>Quality Profile</b>	<b>Specification</b>
Appearance (Color)	Very Dark Red and Very Dark Red-Brown
Appearance (Form)	Powder
Titration with AgNO <sub>3</sub>	≥ 90 %
Purity (HPLC) "Sum of Components Acriflavine 70 %"	5.0 - 95.0 %
Infrared Spectrum	Conforms to Structure
Luminescence	Corresponds to Requirements
Fluorescence (Emmission Max)	485.0 - 495.0 nm

**Intended Use:** For research use only.

This document has been created digitally and is valid without a signature.