



Product Data Sheet

Acryl/Bis solution (19: 1), 40% (w/v)

Product No. GB-CHEM-00658
Volume. 500ML

Quality Profile	Specification
Appearance (Turbidity)	Clear
Appearance (Color)	Colorless
Appearance (Form)	Liquide
Conductivity	≤ 10
Solubility Suitable for use in electrophoresis	Suitable
Heavy Metals (as Lead)	≤ 5 ppm
Iron (Fe) by ICP	≤ 1.0 ppm
Magnesium (Mg) by ICP	≤ 1.0 ppm
Titrateable acid as Acrylic Acid	≤ 0.005 %
Ratio Acrylamide:Bis-Acrylamide ratio is 18:1 to 20:1	Conforms

Intended Use: For research use only.

This document has been created digitally and is valid without a signature.