



Product Data Sheet

Adenine

Product No. GB-CHEM-00659
Volume. 25G

Quality Profile	Specification
Appearance (Color)	White to Off-White
Appearance (Form)	Powder
Solubility (Colour)	Colorless to Faint Yellow
Solubility (Turbidity) 20 mg/mL, 0.5 M HCl	Clear to Slightly Hazy
EMM pH 1.0	12.5 - 13.1
Wavelength	262 - 263 NM
Purity (HPLC)	≥ 99 %

Intended Use: For research use only.

This document has been created digitally and is valid without a signature.