



Product Data Sheet

Adenosine

Product No. GB-CHEM-00660
Volume. 5G

Quality Profile	Specification
Appearance (Color)	White
Appearance (Form)	Powder
Solubility (Colour)	Colorless
Solubility (Turbidity) 50 mg/mL, 1 M NH ₄ OH	Clear
UV/Vis absorbance EMM at pH 7.0	14.6 - 15.2
Wavelength	259 - 260 NM
Purity (HPLC)	≥ 99 %

Intended Use: For research use only.

This document has been created digitally and is valid without a signature.