



Product Data Sheet

Alcian blue 8GX

Product No. GB-CHEM-00668
Volume. 5G

Quality Profile	Specification
Appearance (Color)	Blue to Very Dark Blue and Very Dark Purple Purple to Very Dark Purple
Appearance (Form)	Powder or Crystals
Solubility (Color)	Blue to Very Dark Blue
Solubility (Turbidity)	Conforms to Requirements
Wavelength (1)(UV)	605 - 625 nm
Extinction (1)(UV) (1%/1cm)	≥ 60
Solvent (UV)	Water

Intended Use: For research use only.

This document has been created digitally and is valid without a signature.