



Product Data Sheet

Avidin–Alkaline Phosphatase

Product No. GB-CHEM-00670
Volume. 5KU

Quality Profile	Specification
Appearance (Form)	Liquid
Appearance (Turbidity)	Clear
Appearance (Color)	Almost Colorless to Light Yellow
ELISA - Titer	1:70,000-1:100,000
WB-Cell Line/Tissue Extract	Hela Cells
WB-Western Blot Titer	1:150,000-1:300,000

Intended Use: For research use only.

This document has been created digitally and is valid without a signature.