



Product Data Sheet

Amido Black Stain Solution Ready-To-Use 5X Strength

Product No. GB-CHEM-00671
Volume. 100ML

Quality Profile	Specification
Appearance (Color)	Blue to Black
Appearance (Form)	Liquid
Appearance (Turbidity)	Clear
Suitability Suitable for staining protein bands on nitrocellulose membranes (western blots)	Suitable

Intended Use: For research use only.

This document has been created digitally and is valid without a signature.

Goslar Biotech, 255A Barking Road East Ham, London E6 1LB, United Kingdom
Email: info@goslarbiotech.com, Website: www.goslarbiotech.com