

Product Data Sheet

Ammonium acetate

Product No. GB-CHEM-00673
Volume. 500G

Quality Profile	Specification
Appearance (Color)	White
Appearance (Form) Moist Crystals	Conforms to Requirements
Titration Titration by HCl	≥ 97.0 %
PH at 25°C (c=5%, H ₂ O)	6.7 - 7.3
Residue on Ignition	≤ 0.01 %
Infrared Spectrum	Conforms to Structure
Meets ACS Requirements	Corresponds to Requirements
ACS Specifications Insoluble Matter (10%, H ₂ O) ≤ 0.005 %, Heavy Metals (as Lead) ≤ 5 ppm	Corresponds to Requirements
Iron (Fe)	≤ 5 ppm
Total Sulfur as SO ₄ (ICP)	≤ 0.001 %
Chloride (Cl)	≤ 5 ppm
Nitrate (NO ₃)	≤ 0.001 %

Intended Use: For research use only.

This document has been created digitally and is valid without a signature.