



Product Data Sheet

Amoxicillin trihydrate

Product No. GB-CHEM-00674
Volume. 5G

Quality Profile	Specification
Appearance (Color)	White to Light Yellow to Light Beige
Appearance (Form)	Powder or Crystals
Purity (HPLC area %)	≥ 95.0 %
Specific Rotation (20/D)	248.0 - 268.0°C
Concentration	C=1 In 0.1m Hydrochloric Acid
Melting Point	193 - 203 °C
Water	≤ 15.0 %
Proton NMR Spectrum	Conforms to Structure

Intended Use: For research use only.

This document has been created digitally and is valid without a signature.