



Product Data Sheet

Ampicillin, sodium salt

Product No. GB-CHEM-00675
Volume. 25G

Quality Profile	Specification
Appearance (Form)	Powder
Appearance (Color)	White to Off White
Infrared Spectrum	Conforms to Structure
Water	≤ 2.0 %
Assay (calculated on the anhydrous basis) Endotoxin Level ≤ 0.15 EU/mg	845 - 988 µg/mg
Cell Culture Test	Pass

Intended Use: For research use only.

This document has been created digitally and is valid without a signature.