



Product Data Sheet

AMPSO sodium salt

Product No. GB-CHEM-00676
Volume. 25G

| Quality Profile | Specification |
|---|------------------------|
| Appearance (Color) | White |
| Appearance (Form) | Crystalline Powder |
| Solubility (Color) | Colorless |
| Solubility (Turbidity) 10 g plus 30 ml of H ₂ O | Clear to Slightly Hazy |
| Water (by Karl Fischer) | ≤ 7 % |
| ¹ H NMR Spectrum A290 UV absorbance | Conforms to Structure |
| Heavy Metals | ≤ 5 ppm |
| Titration with H ₂ SO ₄ (Anhydrous) | ≥ 98 % |

Intended Use: For research use only.

This document has been created digitally and is valid without a signature.