

## Product Data Sheet

### Anisomycin

**Product No.** GB-CHEM-00677  
**Volume.** 5MG

Quality Profile	Specification
Appearance (Form)	Powder
Appearance (Color)	White to Faint Yellow
Solubility (Solvent)	Methanol
Solubility (Conc)	19.60 - 20.40 mg/mL
Solubility (Turbidity)	Clear
Solubility (Color)	Colorless to Faint Yellow
NMR (Solvent)	CD3OD
Identity by NMR	Consistent
Purity (HPLC)	≥ 98.00 %
Suitability (Method)	Solubility
Suitability (Conc.mg/ml)	90.00 - 120.00 mg/mL
Suitability (Final Solvent)	Saline, PH 7
Suitability (Clarity)	Clear
Suitability (Color)	Colorless to Faint Yellow

**Intended Use:** For research use only.

This document has been created digitally and is valid without a signature.