



Product Data Sheet

D-Arabinose

Product No. GB-CHEM-00680
Volume. 25G

Quality Profile	Specification
Appearance (Color)	White to Off-White
Appearance (Form)	Powder
Solubility (Color)	Colorless to Faint Yellow
Solubility (Turbidity) 50mg/mL, H ₂ O	Clear
Specific Rotation C = 4 in H ₂ O at 20° C	-105 - -103 °
Purity (GC)	≥ 98 %

Intended Use: For research use only.

This document has been created digitally and is valid without a signature.