



Product Data Sheet

Bacitracin zinc

Product No. GB-CHEM-00684
Volume. 5G

Quality Profile	Specification
Appearance (Form)	Powder
Appearance (Color)	Off White
Zinc (Zn)	4.0 - 6.0 %
Loss on Drying	≤ 5.0 %
Identification Meets the requirements of the test for Composition of Bacitracin	Conforms
Assay (on the dried basis) Bacitracin Units/mg	≥ 65

Intended Use: For research use only.

This document has been created digitally and is valid without a signature.