



Product Data Sheet

Barium acetate

Product No. GB-CHEM-00686

Volume. 250G

Quality Profile	Specification
Appearance (Color)	White
Appearance (Form)	Powder or Crystal or Chunk(s)
ICP Major Analysis Confirms Barium Component	Confirmed
Purity 99.999% Based On Trace Metals Analysis	Meets Requirements
Trace Metal Analysis	≤ 15.0 ppm
Calcium (Ca)	2.0 ppm
Iron (Fe)	0.8 ppm
Lithium (Li)	0.5 ppm
Lead (Pb)	0.5 ppm
Strontium (Sr)	1.0 ppm
Zinc (Zn)	0.5 ppm

Intended Use: For research use only.

This document has been created digitally and is valid without a signature.