



Product Data Sheet

Bromocresol green free acid

Product No. GB-CHEM-00693
Volume. 5G

Quality Profile	Specification
Appearance (Form)	Powder
Infrared Spectrum	Conforms to Structure
Titration by AgNO ₃ After Oxygen Combustion	≥ 95.0 %
Clarity of Solution	Pass
Visual Interval pH 3.8 Yellow	Pass
Visual Interval pH 5.4 Blue	Pass

Intended Use: For research use only.

This document has been created digitally and is valid without a signature.