



Product Data Sheet

Calcium acetate monohydrate

Product No. GB-CHEM-00695
Volume. 250G

Quality Profile	Specification
Appearance (Color)	White
Appearance (Form)	Powder
ICP Major Analysis Confirms Calcium Component	Confirmed
Complexometric EDTA	≥ 99.0 %
Alkalinity	Pass
Insoluble Matter c = 10%; Water	≤ 0.005 %
Barium	≤ 0.01 %
Chloride Content	≤ 0.001 %
Sulfate (SO ₄)	≤ 0.01 %
Magnesium (Mg)	≤ 0.05 %
Iron (Fe)	≤ 0.001 %
Sodium (Na)	≤ 0.02 %
Potassium (K)	≤ 0.01 %
Strontium (Sr)	≤ 0.05 %
Titrateable Acid (meq/g)	≤ 0.035 %

Intended Use: For research use only.

This document has been created digitally and is valid without a signature.