



## Product Data Sheet

### CHES

**Product No.** GB-CHEM-00697  
**Volume.** 500G

Quality Profile	Specification
Appearance (Color)	White
Appearance (Form)	Crystalline Powder
Solubility (Color)	Colorless
Solubility (Turbidity) at 15 g plus 85 ml of H <sub>2</sub> O	Clear
Infrared Spectrum	Conforms to Structure
Water (by Karl Fischer)	≤ 1 %
A <sub>290</sub> UV absorbance 15% W/W	≤ 0.05
ICP-MS (Total: Ag, As, Bi, Cd, Cu, Hg, Mo, Pb, Sb, Sn)	≤ 5 ppm
Titration with NaOH Anhydrous	≥ 99.0 %

**Intended Use:** For research use only.

This document has been created digitally and is valid without a signature.