



Product Data Sheet

Chloramphenicol

Product No. GB-CHEM-00698
Volume. 50G

Quality Profile	Specification
Appearance (Color)	White to Faint Yellow to Faint Brown
Appearance (Form)	Powder or Crystals
PURITY (HPLC AREA %)	≥ 98.0 %
Specific Rotation (20/D)	18.0 - 22.0 DEGREES
Concentration	C=4.8 IN ETHANOL
Melting Point	148 – 154°C
Water	≤ 1.0 %
Proton NMR Spectrum	Conforms to Structure

Intended Use: For research use only.

This document has been created digitally and is valid without a signature.