



## Product Data Sheet

### Chloroform (Trichloromethane)

**Product No.** GB-CHEM-00701  
**Volume.** 1L

Quality Profile	Specification
Appearance (Clarity)	Clear
Appearance (Color)	Colorless
Appearance (Form)	Liquid
Proton NMR Spectrum	Conforms to Structure
Purity (GC)	≥ 99.50 %
Amylenes as Stabilizer	100 - 200 ppm

**Intended Use:** For research use only.

This document has been created digitally and is valid without a signature.