



Product Data Sheet

Citric acid, monohydrate

Product No. GB-CHEM-00704
Volume. 500G

Quality Profile	Specification
Appearance (Color)	White
Appearance (Form)	Powder
Solubility (Color)	Colorless
Solubility (Turbidity) 100 mg/mL, H ₂ O	Clear
Water (by Karl Fischer)	≤ 10.5 %
¹ H NMR Spectrum	Conforms to Structure
Titration with NaOH	≥ 98 %
Purity (GC)	≥ 98 %

Intended Use: For research use only.

This document has been created digitally and is valid without a signature.