

Product Data Sheet

Colchicine

Product No. GB-CHEM-00705
Volume. 100MG

Quality Profile	Specification
Appearance (Color)	White to Yellow
Appearance (Form)	Powder
Solubility (Color)	Yellow to Yellow-Green
Solubility (Turbidity) 50 mg/mL, EtOH	Clear to Slightly Hazy
Solvent Content (by NMR) Ethyl Acetate	≤ 10 %
¹ H NMR Spectrum	Conforms to Structure
¹³ C NMR Spectrum	Conforms to Structure
Purity (HPLC)	≥ 95 %

Intended Use: For research use only.

This document has been created digitally and is valid without a signature.