



Product Data Sheet

o-Cresolphthalein complexon

Product No. GB-CHEM-00709
Volume. 25G

Quality Profile	Specification
Appearance (Color)	White to Off White
Appearance (Form)	Powder
Solubility (Color)	Dark Purple
Solubility (Turbidity) 1 mg/ml, MeOH	Clear
Infrared Spectrum	Conforms to Structure
Water (by Karl Fischer)	≤ 10 %
Carbon (anhydrous)	58.0 - 62.8 %
Nitrogen (anhydrous)	4.0 - 4.7 %
Suitability Suitable for use in Calcium determination	Suitable

Intended Use: For research use only.

This document has been created digitally and is valid without a signature.

Goslar Biotech, 255A Barking Road East Ham, London E6 1LB, United Kingdom
Email: info@goslarbiotech.com, Website: www.goslarbiotech.com