

## Product Data Sheet

### **AGAR AGAR**

**Product No. GB-DCM-00002-1A**

#### **Product Description:-**

<b>Test</b>	<b>Specification</b>
Appearance (Colour) Light Yellow to Beige	Conforms
Appearance (Form)	Powder
Residue on Ignition (Ash)	2.0 - 4.5 %
Loss on Drying	≤ 20 %
Congeaing Temperature	≤ 38 °C
1.5%, H <sub>2</sub> O	
pH	6.0 - 8.0
1.5%, H <sub>2</sub> O, room temperature	
Heavy Metals	Pass
< / = 0.004%	
Cell Culture Testing	Pass
Plant Cell Culture Test	Pass

#### **Precautions and Disclaimer**

This product is for research use only.