



Product Data Sheet

FDA BROTH (ATCC BACTERIOSTASIS BROTH)

Product No. GB-DCM-00213-1A

INTENDED USE

For routine antimicrobial testing of antiseptics and disinfectants.

PRODUCT SUMMARY

ATCC Bacteriostasis Broth (FDA Broth) is useful for subcultures in phenol coefficient and dilution tests of bacteriostatic, germicidal, sporicidal activity and also as a base for the preparation of ATCC Bacteriostasis Agar.

Product Specifications

Ingredients	Gms / Ltr
Peptic digest of animal tissue	10.000
Beef extract	5.000
Sodium chloride	5.000

PRINCIPLE

The medium consists of Peptic digest of animal tissue and beef extract that are sources of carbon, nitrogen, vitamins and minerals. Sodium chloride provides essential ions.

INSTRUCTION FOR USE

- Dissolve 20.0 grams in 1000 ml purified / distilled water.
- Heat if necessary to dissolve the medium completely.
- Sterilize by autoclaving at 15 psi pressure (121°C) for 15 minutes.
- Mix well and dispense as desired.

QUALITY CONTROL SPECIFICATIONS

Appearance of Powder: Cream to yellow homogeneous free flowing powder.

Appearance of prepared medium: Amber coloured clear solution in tubes.

pH (at 25°C) : 6.8±0.2



GOSLAR BIOTECH
Your Partner in Sera

Microorganism	ATCC	Inoculum (CFU/ml)	Growth	Incubation Temperature	Incubation Period
Escherichia coli	25922	50-100	Good-luxuriant	35-37°C	18-24 Hours
Staphylococcus aureus subsp. aureus	25923	50-100	Good-luxuriant	35-37°C	18-24 Hours
Staphylococcus aureus	6538	50-100	Good-luxuriant	35-37°C	18-24 Hours
Salmonella Typhi	6539	50-100	Good-luxuriant	35-37°C	18-24 Hours

STORAGE

Dehydrated powder, hygroscopic in nature, store in a dry place, in tightly-sealed containers between 25-30°C and protect from direct sunlight. Under optimal conditions, the medium has a shelf life of 4 years. When the container is opened for the first time, note the time and date on the label space provided on the container. After the desired amount of medium has been taken out replace the cap tightly to protect from hydration.

Product Deterioration: Do not use if they show evidence of microbial contamination, discoloration, drying or any other signs of deterioration.

DISPOSAL

After use, prepared plates, specimen/sample containers and other contaminated materials must be sterilized before discarding.

This product is for research use only.