

Product Data Sheet

SCHAEDLER AGAR

Product No. GB-DCM-00457-1A

INTENDED USE

For enumeration of various aerobic and anaerobic bacterial species in gastrointestinal tract.

PRODUCT SUMMARY

Schaedler Agar was originally formulated by Schaedler et al and further modified by Mata et al with formulation changes for cultivation and enumeration of aerobic and anaerobic microorganisms. Schaedler Agar supplemented with Vitamin K1 and 5% sheep blood is used for the recovery of fastidious anaerobic bacteria such as Bacteroides. Inclusion of Colistin and Nalidixic acid in the formulation (Schaedler CNA Agar) along with 5% sheep blood is used for the selective isolation of the anaerobic gram-positive cocci, Peptococcus and Peptostreptococcus species. Inclusion of Kanamycin and Vancomycin in the formulation (Schaedler KV Agar) along with 5% sheep blood is used for selective isolation of gramnegative anaerobes. Schaedler Agar serve as an excellent basal media to which blood or other enrichments can be added to enhance the recovery of fastidious anaerobic organisms.

Product Specifications

Ingredients	Gms / Ltr
Tryptone	5.670
Proteose peptone	5.000
Soya peptone	1.000
Yeast extract	5.000
Dextrose (Glucose)	5.830
Sodium chloride	1.670
Dipotassium hydrogen phosphate	0.830
Tris (hydroxymethyl) aminomethane	3.000
L-Cystine	0.400
Hemin	0.010
Agar	15.000

PRINCIPLE

The medium consists of Tryptone, Proteose peptone and Soya peptone, Yeast extract and L-cystine which provide nitrogenous growth factors, vitamins and other essential growth nutrients. Dextrose serves as energy source. Hemin and sheep blood stimulates the growth of fastidious microorganisms and stimulates growth of other Bacteroides species and gram-positive spore formers. Addition of Sodium Polyanethol Sulphonate (SPS) is recommended when using this medium for blood culture. It inhibits phagocytosis and neutralizes the antibacterial activity of fresh blood components. Vitamin K1 enables the cultivation of Bacteroides melaninogenicus and stimulates growth of other Bacteroides species and grampositive spore formers.

INSTRUCTION FOR USE

- Dissolve 43.41 grams in 950 ml purified/distilled water.
- Heat to boiling with frequent agitation to dissolve the medium completely.
- Sterilize by autoclaving at 15 psi pressure (121°C) for 15 minutes.
- Cool to 45-50°C and add 5% sterile defibrinated blood if desired.
- Mix well before dispensing. Avoid overheating and photo-oxidation of the medium, as it will retard the growth of bacteria.
- If desired, add rehydrated contents of 1 vial each of Vitamin K1 Supplement and CNA Supplement to prepare Schaedler CNA Agar or to prepare Schaedler KV Agar, aseptically add rehydrated contents of 1 vial each of Vitamin K1 Supplement and KV Supplement respectively to 1000 ml of Schaedler Agar.

Microorganism	ATCC	Inoculum (CFU/ml)	Growth	Recovery	Incubation Temperature	Incubation Period
Bacteroides fragilis	25285	10-100	luxuriant	20-30%	35-37°C	18-48 Hours
Clostridium sporogenes	13732	50-100	luxuriant	20-30%	35-37°C	18-48 Hours
Clostridium perfringens	12924	50-100	luxuriant	>=50%	35-37°C	18-48 Hours
Clostridium sporogenes	11437	50-100	luxuriant	>=50%	35-37°C	18-48 Hours
Escherichia coli	25922	50-100	Inhibited	0-10%	35-37°C	18-48 Hours
Streptococcus pyogenes	19615	50-100	luxuriant	20-40%	35-37°C	18-48 Hours



QUALITY CONTROL SPECIFICATIONS

Appearance of Powder : Cream to yellow homogeneous free flowing powder.

Appearance of prepared medium: Light amber coloured clear to slightly opalescent gel forms in Petri plates.

pH (at 25°C) : 6.3± 0.2

STORAGE

Dehydrated powder, hygroscopic in nature, store in a dry place, in tightly-sealed containers between 25-30°C and protect from direct sunlight. Under optimal conditions, the medium has a shelf life of 4 years. When the container is opened for the first time, note the time and date on the label space provided on the container. After the desired amount of medium has been taken out replace the cap tightly to protect from hydration.

Product Deterioration: Do not use if they show evidence of microbial contamination, discoloration, drying or any other signs of deterioration.

DISPOSAL

After use, prepared plates, specimen/sample containers and other contaminated materials must be sterilized before discarding.

This product is for research use only.