



## Product Data Sheet

**Horse Serum, Heat Inactivated, Sterile-Filtered**  
**Product No. GB-HS-500A**  
**Packing 500 ml**

### **Product Description**

Horse Serum contains more immunoglobulins than Fetal Bovine Serum and has increased protein content. They are a cost-effective alternative to Fetal Bovine Serum. Horse Serum features a controlled collection of the raw serum. In contrast to slaughtered horses, horse serum has lower endotoxin content and exhibits a lower level of batch to batch variation.

### **Features:-**

- ✚ Alternative to Calf Serum
- ✚ Suitable for ES cell culture
- ✚ Autologous cell culture system
- ✚ Free from fungal and other contaminations

### **Product Specifications**

Appearance	Clear amber liquid
Storage and shelf life	Store at $\leq -15^{\circ}\text{C}$ . Avoid repeated freeze-thaw cycles. Preparation of aliquots recommended. Once opened store at $4^{\circ}\text{C}$ and use within 4-6 weeks.
Shipping conditions	Frozen (Dry ice)
Thawing	Overnight at $+2^{\circ}\text{C}$ to $+8^{\circ}\text{C}$ . Swirl gently to homogenize

### **Precautions and Disclaimer**

This product is for research use only.