



GOSLAR BIOTECH
Your Partner in Sera

Product Data Sheet

(Other Animal Sera, Heat Inactivated)

Product No. GB-PS-500A, GB-PS-100A – Porcine Serum

Product No. GB-SS-500A, GB-SS-100A – Sheep Serum

Product No. GB-GS-500A, GB-GS-100A – Goat Serum

Product No. GB-MS-100A – Mouse Serum

Product No. GB-RS-100A – Rabbit Serum

Product No. GB-RAT-100A – Rat Serum

Product Description

Other animal sera are used in various applications, e.g. autologous cell culture systems or immune-histological methods.

Porcine serum is often used in cultures that proliferate mycoplasma or viruses. Goat and Sheep Serum are used in immunohistochemistry and to eliminate non-specific binding reactions.

Applications

- Autologous cell culture systems
- Immunohistochemistry
- Immunoblotting
- Immunofluorescence
- Immunoprecipitation

Test Parameter

- pH
- Osmolality
- Endotoxin
- Total protein
- Hemoglobin
- Mycoplasma
- Sterility Storage: Store at $\leq -20^{\circ}\text{C}$

Precautions and Disclaimer

This product is for research use only.